

AD-A069 500

PRINCETON UNIV N J DEPT OF CHEMICAL ENGINEERING  
PYROLYSIS-MOLECULAR WEIGHT CHROMATOGRAPHY-VAPOR PHASE INFRARED --ETC(U)  
MAY 79 E KIRAN, J K GILLHAM  
N00014-76-C-0200  
NL

UNCLASSIFIED

| OF |

AD  
A069500



AD A 069500

DDC FILE COPY

OFFICE OF NAVAL RESEARCH

Contract N00014-76-C-0200

Task No. NR 356-504

TECHNICAL REPORT NO. 16

12  
6.5.

Pyrolysis-Molecular Weight Chromatography-Vapor  
Phase Infrared Spectrophotometry: An On-Line  
System for Analysis of Polymers

by

E. Kiran and J. K. Gillham

Prepared for Publication as a Chapter in


"Developments in Polymer Degradation-2"

Editor: Professor N. Grassie

Applied Science Publishers Ltd., London, England

Polymer Materials Program  
Department of Chemical Engineering  
Princeton University  
Princeton, NJ 08544

May 1979



Reproduction in whole or in part is permitted for  
any purpose of the United States Government

This document has been approved for public release  
and sale; its distribution is unlimited

J. K. Gillham  
Principal Investigator  
609/452-4694

79 06 04 103

## PYROLYSIS-MOLECULAR WEIGHT CHROMATOGRAPHY-VAPOR PHASE

## INFRARED SPECTROPHOTOMETRY: AN ON-LINE SYSTEM

## FOR ANALYSIS OF POLYMERS

Erdoğan Kiran  
SEKA (Türkiye Selüloz ve Kağıt Fabrikaları İşletmesi)  
Central Research Laboratory  
İzmit, Kocaeli, Turkey

and

John K. Gillham  
Polymer Materials Program  
Department of Chemical Engineering  
Princeton University  
Princeton, New Jersey 08540 USA

Accession For	
NTIS GRA&I	<input checked="" type="checkbox"/>
DDC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By _____	
Distribution/	
Availability Codes	
Dist	Avail and/or special
A	

An instrumental system consisting of a combination of a programmable pyrolyzer, a thermal conductivity detector, a mass chromatograph, a gas chromatograph, a fast-scan vapor-phase infrared spectrophotometer and a computer is discussed and reviewed with examples which show the utility of the system in the analysis of decomposition phenomena in polyolefins, polyolefinsulfones, polymethacrylates, polystyrenes, and polybutadienes.

The system records the thermal history before and during pyrolysis of the sample and provides chromatographic retention times, infrared spectra, mass numbers and relative amounts of the volatile products of pyrolysis. Identification of the constituents of the effluents is thus facilitated and mechanisms for decomposition become easier to propose and to verify.

## INTRODUCTION

Pyrolytic decomposition followed by analysis of the volatile products of pyrolysis is a powerful method for characterization of polymeric materials. The scope and applicability of this approach are determined by the type of pyrolyzer and the separation and identification techniques that are utilized.

The on-line system developed at Princeton University comprises a programmable pyrolyzer, a thermal conductivity cell, a trap, a mass chromatograph, a conventional gas chromatograph, a vapor-phase infrared spectrophotometer, and a digital computer (1,2,3). The system has proved powerful for analysis of polymers as has been shown by results obtained with polyolefins (1,4,5,6,7), polyolefinsulfones (1,8), polymethacrylates (1), polystyrenes (1,3,8,9) and polybutadienes (10). This report provides a brief description of the system and reviews results obtained with these polymers.

## THE PYROLYSIS LABORATORY

A schematic diagram of the system is shown in Figure 1. A critical discussion of its various parts has been presented elsewhere (2). Briefly, the system consists of a programmable pyrolyzer, a thermal conductivity cell (TC), a trap, a mass chromatograph (MC), a conventional gas chromatograph (GC), a vapor-phase infrared spectrophotometer (Vapor Phase IR) and a digital computer (IBM System 7).

Figure 2 shows details of the coupling of the pyrolyzer with the thermal conductivity detector, the trap, and the mass chromatograph.

The pyrolyzer consists of a quartz tube and an external heater which provides flexibility for either program-heating or flash-heating.



The thermal conductivity cell monitors the formation of the volatile products of pyrolysis as a function of temperature and time. It provides immediate information about the onset and the progress of the thermal events that result in the formation of volatile products. The response curves contain features that are characteristic of the material under investigation.

The combination of the pyrolyzer with the thermal conductivity detector is complementary to a thermogravimetric analyzer and/or differential thermal analyzer. As an example, the thermal history before and during pyrolysis of poly(pentene-1-sulfone) as recorded with the thermal conductivity cell (TC), a thermogravimetric analyzer (TGA) and a differential thermal analyzer (DTA) is shown in Figure 3.

The trap is an important part of the interface between the pyrolysis assembly and the chromatographic assembly. It is a short stainless steel column packed with Porapak Q (Waters Associates) and has a geometry which permits a fraction of it to be placed in a Dewar flask for subambient cooling. The function of the trap is to collect the fragments of pyrolysis that are formed in the temperature range of interest (i.e. selective trapping). The products are released from the trap by rapid heating, and introduced as a slug to the chromatographic columns for separation. The valving system permits monitoring of the release of the products from the trap with the same thermal conductivity cell used to monitor the formation of volatile products; the response is recorded as a sharp peak, as shown for example in Figure 4.

The mass chromatograph is a special type of a gas chromatograph which directly provides mass numbers of the resolved components of a mixture by

means of a pair of gas density detectors (11). The instrument consists of two independent gas chromatographs which use a common injection port. When a sample is introduced to the instrument it is split into two approximately equal fractions (splitter, Figure 2) and each fraction is carried to a trap which is similar to the trap associated with the pyrolyzer. Then, through a valve arrangement, different carrier gases are introduced to the two traps; the samples, released by rapid heating of the traps, are carried over to two matched chromatographic columns for separation. The constituents eluting from the chromatographic columns are detected by gas density detectors. The recorder output from the mass chromatograph displays two sets of peaks corresponding to the responses from the two gas chromatographic systems with different carrier gases for the same constituents of the mixture (see Figures 3 and 4).

Theoretical aspects and operational principles of the mass chromatograph and the gas density detector are presented elsewhere (1,11). The ratio of the responses is related to the molecular weight ( $M_x$ ) of a constituent through

$$M_x = \frac{\left(\frac{A_1}{A_2}\right) K M_{C_2} - M_{C_1}}{\left(\frac{A_1}{A_2}\right) K - 1}$$

where  $M_{C_1}$  and  $M_{C_2}$  are the molecular weights of the carrier gases 1 and 2,  $A_1$  and  $A_2$  are the chromatographic peak areas for the constituent, and  $K$  is the instrument constant determined from analysis of samples of known molecular weights from

$$K = \left( \frac{M_x - M_{C_1}}{M_x - M_{C_2}} \right) / \left( \frac{A_1}{A_2} \right).$$

The response of the gas density detector is proportional to  $(M_x - M_c)/M_x$  and therefore, for a given carrier gas of molecular weight  $M_c$ , the response tends to a constant as  $M_x$  increases; the lower the molecular weight of the carrier gas, the faster is this tendency. In order to cover a wide range of molecular weights, the two carrier gases are chosen to be of highly different molecular weights. The gases used in these studies were  $\text{CO}_2$  and  $\text{ClC}_2\text{F}_5$  (Freon-115). Use of  $\text{CO}_2$  and  $\text{SF}_6$  is reported (12) in a study of oxidative degradation of isotactic poly(1-pentene) using a combination of a thermal conductivity detector and a mass chromatograph.

Since the molecular weights of the carrier gases are known and since the response ratios  $(A_1/A_2)$  can be measured from the chromatographic output, calculation of the molecular weights of the constituents becomes a simple process after the instrument constant is evaluated. In principle, the instrument constant  $K$  can be evaluated by injecting just one solute (of known molecular weight) into the mass chromatograph and measuring the response ratio from the output. However, the accuracy of measuring the response ratio varies with molecular weight and consequently  $K$  must be evaluated by using more than one compound if it is to be applicable to a wide range of molecular weights. Otherwise instantaneous  $K$  values applicable in small ranges are used.

In addition to providing mass numbers of resolved constituents the mass chromatograph can be used to provide estimates of the weight amounts of each ( $w_i$ ). This is performed by comparing the areas on the  $\text{CO}_2$  channel output of the mass chromatograph with those of peaks obtained under the same chromatographic conditions from a synthetic mixture containing known amounts of n-hydrocarbons (3). For compounds of very similar molecular



weight and with a similar retention time, responses of the gas density balance detector are directly proportional to their amounts (see equation 71 of reference 11). The relative amount of a component as a weight percentage ( $P_{wt}$ ) of those measured can be calculated using  $P_{wt} = (w_i / \sum w_i) \times 100$ .

The infrared spectrophotometer is a fast-scan vapor-phase infrared spectrophotometer and is suitable for taking "on-the-fly" spectra of constituents eluting from chromatographic columns. The instrument is capable of scanning from 2.5 to 15  $\mu\text{m}$  (4000 to 670  $\text{cm}^{-1}$ ) in either 6 or 30 seconds. The spectrophotometer can be coupled with either the mass chromatograph or the conventional gas chromatograph. When coupled with the mass chromatograph sensitivity is reduced due to the fact that column effluents are highly diluted by the reference gas entering the gas density detector. Furthermore the IR spectra of the carrier gases,  $\text{CO}_2$  and Freon-115, are complex and can mask spectra of constituents (2). Therefore coupling of infrared spectrophotometer with a conventional gas chromatograph is preferred. A schematic of the coupling of the gas chromatograph with the infrared spectrophotometer and the computer is shown in Figure 5.

The gas chromatograph used is a computer-compatible, automated research instrument equipped with thermal conductivity and flame ionization detectors.

The computer facility is a flexible mini-computer capable of pre-processing of data for analysis by a much larger computational system (to which it is connected by telephone). As a constituent elutes from the chromatographic column it may be trapped in the IR cell and a set number of IR spectra taken. These spectra are computer-averaged and the background spectrum is subtracted in order to increase the signal-to-noise ratio.



## THERMAL DECOMPOSITION OF POLYMERS

Thermal decomposition of various polymers has been studied by the on-line system just described. Polymers studied include:

1) polyolefins (1,4,5,6,7); 2) polyolefinsulfones (1,8); 3) polymethacrylates (1); 4) polystyrenes (1,3,8,9) and 5) polybutadienes (10).

Comparative discussions with literature reports on the same polymers using different techniques have been presented in references 1 to 11.

1) Polyolefins



Polyethylene, polypropylene, and polyisobutylene samples were pyrolyzed in flowing helium at 20°C/min heating rate from room temperature to 600°C. Products volatilizing in the temperature range 300-600°C were selected for trapping and analysis.

The output from the mass chromatograph for low density polyethylene is shown in Figure 6. The regularly spaced doublet peaks were shown to correspond to alkenes and alkanes by calculating the molecular weights and comparing the retention times with hydrocarbons (1,7). Molecular weights calculated from the mass chromatographic output for some of the assigned structures are given in Table I. The formation of the hydrocarbon constituents seen in Figure 6 and the observed relative abundance of  $C_6$ ,  $C_{10}$ ,  $C_{14}$  and  $C_{18}$  alkenes and  $C_3$ ,  $C_7$ ,  $C_{11}$  and  $C_{15}$  alkanes were considered to follow from an intramolecular radical transfer process which involves the 5th carbon atom of the primary macroradical through a pseudo six-membered ring intermediate. The presence of ethylene is not seen in Figure 6 since the column conditions were not suitable for its analysis.


Main chain scission in polyethylene results initially in primary radicals which can undergo either depolymerization or intra and/or intermolecular radical transfer processes. The relative extents of the

competing reactions determines the degradation pattern.


Intramolecular radical transfer reactions in primary radicals (I) lead to the formation of secondary radicals (II). These can undergo  $\beta$ -scission to the left of the radical center and produce alkenes (III).



If the  $\beta$ -scission is to the right of the radical center short chain primary radicals (IV) form




which can abstract a hydrogen atom from another molecule and produce alkanes (V).




Highly strained conformations will not be permitted in the intramolecular backbiting steps. The process is therefore argued to involve preferentially the fifth carbon atom through a pseudo six-membered ring intermediate. If subsequent repetitive intramolecular transfer follows before chain scission, the consequence of preferential radical transfer to the 5th, 9th, 13th, 17th, and 21st,..... carbon atoms is the preferential

formation of  $C_6$ ,  $C_{10}$ ,  $C_{14}$ ,  $C_{18}$ ,  $C_{22}$ ,... alkenes (reaction II  $\rightarrow$  III) and  $C_3$ ,  $C_7$ ,  $C_{11}$ ,  $C_{15}$ ,  $C_{19}$ ,... alkanes (reactions II  $\rightarrow$  IV  $\rightarrow$  V) respectively. The probability of sequential migration of the free radical along the chain without undergoing  $\beta$ -scission reactions is expected to decrease as the number of transfer steps increases. One and two step processes should be most dominant and, as observed in Figure 6, among alkenes  $C_6$  and  $C_{10}$ , and among alkanes  $C_3$  and  $C_7$  represent major constituents.



The mass chromatographic output obtained for isotactic polypropylene is shown in Figure 7. Calculation of the molecular weights has shown that the peaks in the chromatogram have molecular weights corresponding to those of the monomer, dimers, trimers and higher oligomers (1, 7). Calculated molecular weights for some of the peaks are included in Table I. The monomer (peak 1), which is formed extensively, is produced from either the primary or the secondary radicals formed on chain homolysis. Intramolecular radical transfers to the 6th, 10th, and 12th carbon atoms in the primary macroradicals and to the 5th and 9th and also the 13th carbon atoms in the secondary macroradicals (indexing from the secondary carbon radical at the chain end) account for the other major products of decomposition. The formation of the other volatile products of pyrolysis of polypropylene and their relative abundances have been explained also in terms of intramolecular radical transfer processes. In polypropylene initiation by chain scission results in primary (I) and secondary (II) radicals.








Intramolecular transfer reactions followed by  $\beta$ -scission to the left of the resulting tertiary radical produces alkenes (III) from primary radicals and alkenes (IV) from secondary radicals.



$\beta$ -scission to the right of the tertiary radical, followed by hydrogen atom abstraction, leads to the formation of alkanes (V) and (VI) from primary and secondary macroradicals, respectively.




Specific major decomposition products, i.e., pentane (peak 3, structure VI,  $n=1$ ), 2,4 dimethyl-1-heptene (peak 9, structure IV,  $n=1$ ) and 2,4,6,8 tetramethyl-1-hendecene (peak 12, structure IV,  $n=3$ ) can form through intramolecular radical transfer in the secondary radicals, and 2-methylpentane




(peak 5, structure V,  $n=2$ ), 2,4,6-trimethylnonane (peak 11, structure V,  $n=4$ ) and 2,4,6,8-tetramethylhendecane (peak 13, structure V,  $n=5$ ) can form from intramolecular radical transfer in the primary radicals. Among the higher molecular weight fragments 2,4,6,8,10,12 hexamethyl-1-tridecene (peak 16) and 2,4,6,8,10,12 hexamethyl-1-pentadecene (peak 17) are formed in larger amounts. Their formation is a continuation of the processes leading to peaks 12 and 13 which are major decomposition products.

The mass chromatogram for polyisobutylene displays a series of regularly spaced groups of peaks which is characteristic of a homologous series of products (Figure 8). In addition to depolymerization processes which account for the extensive formation of the monomer (peak 1), intramolecular radical transfer processes in the primary and tertiary macroradicals which form on initial chain scission account for the formation of dimers, trimers, and high oligomers (1,7). More intramolecular transfer occurs in the primary radicals than the tertiary radicals for energetic reasons. A major product, the trimer 2,4,4,6,6-pentamethyl-1-heptene (peak 12), is formed from a primary radical via intramolecular radical transfer to a methyl group:



When the radical thus formed undergoes  $\beta$ -scission to the left of the adjacent carbon, alkenes with the structure (I) are formed:



For  $n = 2$  the alkene is the trimer 2,4,4,6,6-pentamethyl-1-heptene. In each quartet (peaks 14-29) the first two peaks correspond to alkenes which arise from tertiary radicals, but the last two, which are formed in larger amounts, to alkenes which are formed from primary radicals. For example in structure (I),  $n = 4$  corresponds to peak 20,  $n = 5$  to peak 24. Peaks 21 and 25 arise from primary radicals but via intramolecular radical transfer to methylene instead of a methyl group. Assigned structures and calculated molecular weights for some of the peaks are given in Table 1.

## 2) Polyolefinsulfones

Polybutene-1-sulfone, polypentene-1-sulfone and polyhexene-1-sulfone have been studied by pyrolyzing samples in helium at a heating rate of  $20^\circ\text{C}/\text{min}$  (1,8). The volatile products formed in the temperature range  $100\text{--}500^\circ\text{C}$  were analyzed by the mass chromatograph.

The mass chromatographic outputs are shown for polypentene-1-sulfone in Figure 3. The results for the three polymers are compared in Figure 9 in terms of the outputs of the  $\text{CO}_2$  channel of the mass chromatograph. The primary products of decomposition have been found to be  $\text{SO}_2$  (peak 1) and the respective olefin, i.e., butene, pentene, and hexene (peaks 2). These polysulfones contain  $-\text{C}-\overset{\text{O}}{\underset{\text{O}}{\text{S}}}-$  linkages in the polymer backbone. The bond dissociation energy for the C-S bond (55-60 kcal/mole) is appreciably lower than for the C-C bond (80-85 kcal/mole) and the C-H bond (90-100 kcal/mole) (13). Consequently thermal decomposition in these polymers

may occur readily by scission of the C-S bond and results in the formation of  $\text{SO}_2$  and the respective olefin.

As seen in Figure 9, in going from one chromatogram to the next, the locations of the constituents (peaks 2 to 6) other than  $\text{SO}_2$  (peak 1) shift to higher retention times in an ordered fashion. Calculations of the molecular weights have shown that in each chromatogram the incremental increase in molecular weight in going from peak 2 to 3 is 32 mass units, from peak 3 to 4 is equal to the molecular weight of the respective olefin, from 4 to 5 is 32 mass units, and from 5 to 6 is 16 mass units. Based on the atomic weight of S being 32 and that of O being 16, some compositional possibilities were suggested for peaks 3 to 6. These have the form  $[\text{C}_n^=][\text{S}]$  for peaks 3,  $[\text{C}_n^=]_2[\text{S}]$  for peaks 4,  $[\text{C}_n^=]_2[\text{S}][\text{O}]_2$  or  $[\text{C}_n^=][\text{S}]_2$  for peaks 5, and  $[\text{C}_n^=]_2[\text{S}][\text{O}]_3$  or  $[\text{C}_n^=]_2[\text{S}]_2[\text{O}]$  for peaks 6, where  $\text{C}_n^=$  represents  $\text{C}_4$  olefin,  $\text{C}_5$  olefin or  $\text{C}_6$  olefin. Mechanisms for the formation of these structures were not advanced (18). Their formation probably involves some destruction of  $\text{SO}_2$  in the presence of free radicals.

These polymers all show a two-step decomposition phenomenon (see Fig. 3). Selective trapping and analysis of the volatile decomposition products which form in the temperature range 100–200°C showed that only  $\text{SO}_2$  and the respective olefin are formed in the first step of degradation.

In addition to these polyolefinsulfones, a sample of polystyrene-sulfone has been studied (1,8).  $\text{SO}_2$  and styrene are found to be the major products of decomposition. Additional fragments derived from styrene dimers were observed as indicated by comparison of the results with those obtained from pyrolysis of polystyrenes; this provided evidence for the



presence of consecutive styrene residues in the repeat unit of the polystyrenesulfone.


### 3) Polymethacrylates

Polymethylmethacrylates and polytertiarybutylmethacrylates of different tacticities have been studied by pyrolyzing samples in flowing helium atmosphere at a heating rate of 20°C/min (1).

Thermal conductivity responses during pyrolysis of polymethylmethacrylates showed some differences for the different tactic forms (Figure 10). Analysis of the products of pyrolysis formed in the temperature range 200-500°C shows, however, that all the tactic forms revert almost exclusively to a single compound with a calculated molecular weight which compares well with that of the monomer, methylmethacrylate (Table I). Mass chromatograms are shown in Figure 11. Analysis of the products of pyrolysis formed from the isotactic polymer in the temperature range 420-500°C (see Figure 10) has shown that in this region monomer formation is no longer observed (see Figure 11D). The traces of high molecular weight constituents observed were not characterized.

The thermal histories during pyrolysis of the polytertiarybutylmethacrylates are shown in Figure 12. These polymers behave similarly and show two-stage decomposition. The mass chromatogram of the volatile products of pyrolysis formed in the temperature range 140-500°C is shown in Figure 13 for the isotactic polymer. The major products of decomposition of all of the tactic forms are found to be isobutylene (peak 1) and the monomer, tertiarybutylmethacrylate (peak 2) (Table I). Formation of the monomer is a consequence of depolymerization. The formation of isobutylene is indicative of ester decomposition, i.e.,





The relative extents of these reactions were found to be dependent on the tacticity of the sample as evidenced by differences in the monomer to isobutylene peak height ratios measured in the outputs of the  $\text{CO}_2$  channel of the mass chromatograph (1). Monomer to isobutylene ratio was found to be smallest for the syndiotactic and largest for the isotactic polymer. This observation was interpreted to mean that the extent of depolymerization compared with ester decomposition is greatest with the isotactic polymer; the steric hinderances of the bulky side groups apparently favor depolymerization most in the isotactic polymer. Further analyses of the volatile decomposition products formed in the temperature ranges 140 to 225°C, 140 to 320°C and 320 to 500°C in the isotactic polymer showed that all monomer generation takes place in the temperature range 140 to 320°C. In the second stage of decomposition (i.e., from 320 to 500°C) no monomer is generated, but some isobutylene continues to form.


#### 4) Polystyrenes

Anionic, thermal, and comb-shaped polystyrenes have been studied (3,9).

Anionic polystyrenes with different average molecular weights, which had been polymerized using n-butyl lithium as the initiator and terminated


with protons, were pyrolyzed in flowing helium at a heating rate of 20°C/min, and the products volatilizing in the temperature range 200-500°C were analyzed. The mass chromatographic output for a low molecular weight polymer is shown in Figure 14. Calculations of the molecular weights and examination of the on-the-fly IR spectra showed that the major decomposition products are toluene (peak 4), styrene (peak 5), 1,3-diphenyl-butadiene (peak 11) and 2,4,6 tri-phenyl-1-hexene (peak 16) (Table I and Table II). The relative amounts of these constituents have been found to depend on the molecular weight of the initial polymer. Among the constituents, toluene shows the highest sensitivity to molecular weight (Table II).

In polystyrene cleavage of the backbone results in the primary (I) and secondary (II) radicals.



Both of these can undergo depolymerization and form styrene (peak 5, Tables I and II).


1,3 diphenylbutadiene (peak 11, Tables I and II) forms from tertiary radicals of the form




formed by hydrogen abstraction, which cleave β from the tertiary radical to form terminally unsaturated species with allylic hydrogens:




Abstraction of the allylic hydrogen followed by  $\beta$ -scission of the resulting tertiary radical leads to 1,3 diphenylbutadiene:



Intramolecular radical transfer in secondary radicals of type II, via a pseudo six-membered ring intermediate, may result in tertiary radicals such as:



which upon  $\beta$ -scission to the left of the radical center leads to 2,4,6-triphenyl-1-hexene (peak 16, Tables I and II):



Depolymerization reactions in secondary radical II, if not stopped by other mechanisms, may eventually lead to the primary radical



which can abstract a hydrogen atom from another molecule to yield toluene (peak 4, Tables I and II).



Thermally polymerized polystyrenes have been found to display a different decomposition pattern from that of anionic polystyrenes. Comparison of the mass chromatograms shown in Figures 14 and 15 and examination of Table II shows that fewer products are formed with thermal polystyrenes; those peaks which are highly molecular weight dependent (i.e., influenced by chain ends) are absent.

Another set of low molecular weight anionic polystyrenes initiated with a di-carbanionic species and terminated at both ends with either proton,  $\alpha$ -naphthyl, diphenylmethyl, or hydroxyl groups, was also studied. These polymers have been found to give products of decomposition which are similar to those of thermal polystyrenes. Polymers with  $\alpha$ -naphthyl ends and diphenylmethyl ends showed some differences in that products arising from these groups are also observed. Figure 16 shows the mass chromatogram for the polymer with  $\alpha$ -naphthyl ends. Peaks 5 and 6, predicted to be naphthylethane and 3-naphthyl-1-propene, respectively are characteristic of this polymer.

The major products of pyrolysis of comb-shaped polystyrenes were found to be similar to those of thermal polystyrenes also. However the steric hinderance of the side branches in the comb-shaped polymers was considered to affect the relative abundance of the observed constituents. Formation of the trimer (peak 16) is restricted and lesser amounts form.

#### 5) Polybutadienes


1,4-polybutadienes with different cis/trans ratios have been studied (10).




Thermal histories before and during pyrolysis of a high trans and a high cis 1,4-polybutadiene in an inert atmosphere are shown in Figure 17. The polymers undergo two-step decomposition. Degradation to 15% weight loss and also after 15% weight loss has been investigated. Mass chromatograms of the volatile products formed during 15 weight percent degradation from a high cis and a high trans sample are shown in Figures 18 and 19, respectively. Calculation of the molecular weights and examination of on-the-fly IR spectra (Figure 20) have shown that the major products of decomposition are 1,3-butadiene (monomer, peak 1), cyclopentene (peak 2), 1,3-cyclohexadiene (peak 5) and 4-vinyl-1-cyclohexene (dimer, peak 9). Their relative amounts were found to vary with the cis/trans ratio in the polymer sample. The relative quantities of 4-vinyl-1-cyclohexene and 1,3-butadiene decrease, cyclopentene increases and 1,3-cyclohexadiene is relatively unchanged with increasing trans content as shown in Figure 21.

The mass chromatogram of the products of pyrolysis after 15 percent weight loss of high trans 1,4-polybutadiene (Figure 22) showed that 4-vinyl-1-cyclohexene, which is the most abundant product in degradation to 15% weight loss (see Figure 19), is only a minor product of further decomposition.


The major products of decomposition of polybutadienes are accounted for by a free radical mechanism. Homolytic scission of the single bond of the polymer chain at a site which is  $\beta$  to the double bond results in the radical



The monomer 1,3-butadiene arises from I via  $\beta$ -scission, i.e.,



4-vinyl-1-cyclohexene arises via a backbiting step involving a pseudo six-membered ring intermediate which is favored by a cis 1,4-structure, as is indicated by



1,2-vinyl-polybutadiene has also been studied (10). In contrast to 1,4-polybutadienes, this polymer undergoes one-step decomposition. The mass chromatogram of all of the volatile products of pyrolysis was similar to that from pyrolysis above 15% weight loss for 1,4-polybutadienes (compare Fig. 4 and Fig. 22). This observation has been interpreted to suggest that part of the 1,4-structure undergoes isomerization to the 1,2-structure in the course of being heated to 15% weight loss.

#### COMMENTS ON CALCULATED MASS NUMBERS

The accuracy of the molecular weight determination using a mass chromatograph depends on the invariance of the instrument constant and the accuracy and precision with which the response ratios are measured (11).

Table I shows the theoretical molecular weights and the molecular weights calculated from the mass chromatographic outputs for some of the

assigned structures in the pyrograms discussed in the foregoing sections. As seen from this table, the calculated values are acceptable values, and thus the mass chromatograph offers a convenience not found in the usual gas chromatograph-mass spectrometer combinations. Mass spectrometry, even though more accurate, requires elaborate interfacial systems and sophisticated data-reduction processes.

#### CONCLUDING REMARKS

The results presented in the foregoing sections show the utility and power of the on-line pyrolysis system in polymer analysis. The system provides thermal history before and during pyrolysis and makes selective trapping and analysis of volatile constituents possible. Chromatographic retention times, easily calculable mass numbers, computer-added IR spectra -- all available from an on-line system -- facilitate identification of the volatile constituents. Thus not only are fingerprint pyrograms generated, but also detailed information about the products of pyrolysis and their relative abundances are obtained. Based on these, mechanisms for decomposition become easier to propose and to verify.



## REFERENCES

1. E. Kiran, Ph.D. Thesis, Department of Chemical Engineering, Princeton University, Princeton, N. J. (1974).
2. E. Kiran and J. K. Gillham, J. Appl. Polym. Sci., 20, 931 (1976).
3. H. H. Kuo, Ph.D. Thesis, Department of Chemical Engineering, Princeton University, Princeton, N. J. (1976).
4. E. Kiran and J. K. Gillham, ACS, Polymer Preprints, 14(1), 580 (1973).
5. E. Kiran and J. K. Gillham, Soc. Plastics Engineers, Tech. Papers, 19, 502 (1973).
6. E. Kiran and J. K. Gillham, J. Macromol. Sci.-Chem. A8(1), 211 (1974).
7. E. Kiran and J. K. Gillham, J. Appl. Poly. Sci., 20, 2045 (1976).
8. E. Kiran, J. K. Gillham and E. Gipstein, J. Appl. Poly. Sci., 21, 1159 (1977).
9. H. H. Kuo, H. A. Pfeffer and J. K. Gillham, Amer. Chem. Soc., Coatings and Plastics Preprints, 35(1), 434 (1975).
10. S. Tamura and J. K. Gillham, J. Appl. Poly. Sci., 22, 1867 (1978).
11. E. Kiran and J. K. Gillham, Anal. Chem. 47(7), 983 (1975).
12. S. S. Stivala and S. M. Gabbay, in "Aspects of Degradation and Stabilization of Polymers", H. H. G. Jellinek (ed.), Elsevier, New York, 1978, Chapter 13.
13. T. L. Cottrell, "The Strength of Chemical Bonds", Butterworth, London, 1958.

## ACKNOWLEDGMENT

Partial support of the Office of Naval Research is acknowledged.

TABLE I. Calculated Molecular Weights and Assigned Structures for some of the Characteristic Peaks in the Pyrograms discussed in this paper.

Figure No	Peak No	Assigned Structure	Molecular Weight	
			Theory	Calculated
6	6	Hexene	84	85.1
	6'	Hexane	86	85.5
	10	Decene	140	141.6
	10'	Decane	142	142.9
	14	Tetradecene	196	195.4
	20	Eicosene +(eicosane)	(283)	287.4
7	3	pentane	72	77
	5	2 methylpentane	86	86
	9	2,4 dimethyl-1-heptene	126	124.7
	11	2,4,6 trimethylnonane	170	171
	12	2,4,6,8 tetramethyl-1-hendecene	210	210.7
	16	2,4,6,8,10,12 hexamethyl-1-tridecene	266	270.4
8	6	2,2,4 trimethylpentane	114	112.2
	12	2,4,4,6,6 pentamethyl-1-heptene	168	170.5
	16	2,4,4,6,6,8,8 heptamethyl-1-nonene	224	227
	20	2,4,4,6,6,8,8,10,10 nonamethyl-1-hendecene	280	277.2
	24	2,4,4,6,6,8,8,10,10,12,12 undecamethyl-1-tridecene	336	325
11	1	methylmethacrylate	100	100.3
13	1	isobutylene	56	60.1
	2	tertiarybutylmethacrylate	142	136
14,15	4	toluene	92	93
	5	styrene	104	109
	11	1,3 diphenylbutadiene	206	205
	16	2,4,6 triphenyl-1-hexene	313	339
18,19	1	1,3 butadiene	54	56.3
	2	cyclopentene	68	68.7
	5	1,3-cyclohexadiene	80	77.2
	9	4-vinyl-1-cyclohexene	108	108.8

TABLE II. Volatile Products of Pyrolysis of Anionic and Thermal Polystyrenes. Variation of their Relative Amounts (as weight percentages) with the Molecular Weight of the Initial Polymer.

Peak/ $\bar{M}_n$ Number	Anionic					Thermal		
	2100	4000	10,000.	20,400	390,000	1880	6270	111,000
1	.076	.031	Trace	Trace	0	0	0	0
2	.69	.37	Trace	Trace	0	0	0	0
3	.71	.38	1.77	2.97	1.82	0	0.32	0
4	8.33	5.24	2.79	1.99	1.31	1.0	1.40	2.54
5	37.4	35.8	41.2	45.7	49.0	42.80	48.00	50.90
6	4.14	1.61	Trace	Trace	0	0	0	0
7	2.72	3.64	Trace	Trace	0	0	0	0
8	2.56	1.14	Trace	Trace	0	0	0	0
9	3.02	2.75	3.39	3.04	2.22	0	0	0
10	4.11	5.25	5.75	4.36	5.91	Trace	Trace	Trace
11	8.96	9.76	12.56	11.74	11.37	17.50	18.00	18.87
12	1.66	1.17	.50	Trace	0	0	0	0
13	2.99	2.23	1.27	Trace	0	0	0	0
14	4.53	5.85	6.45	7.19	5.07	0	0	0
15	5.75	7.63	7.38	6.57	8.69	0	0	0
16	22.33	18.14	16.70	16.44	14.64	35.60	29.80	26.33



## FIGURE CAPTIONS

Figure 1. Schematic diagram of the laboratory for pyrolytic studies.

Figure 2. Coupling of the pyrolyzer with the thermal conductivity detector, the trap, and the mass chromatograph.

Figure 3. Left: Thermal history before and during pyrolysis of poly(pentene-1-sulfone) ( $[\eta] = 1.49 \text{ dl/g}$ ).  
Right: Mass chromatogram of the volatile pyrolysis products of poly(pentene-1-sulfone). Peak attenuations were X256 for peaks 1 and 2, and X16 thereafter. The columns (SE-30) were program-heated at  $5^\circ\text{C/min}$  from  $30^\circ$  to  $300^\circ\text{C}$ .

Figure 4. Mass chromatogram of the volatile pyrolysis products of syndiotactic 1,2-polybutadiene (1,2 structure 83-90%,  $\bar{M}_w = 177,000$ ,  $\bar{M}_n = 107,000$ ) formed during 100% weight loss. Peak attenuations for the Freon-115 channel were X64 for peak 1, and X16 thereafter. Those for the  $\text{CO}_2$  channel were X8 for all peaks. The columns (Dexsil 300) were program-heated at  $4^\circ\text{C/min}$  from  $30^\circ$  to  $300^\circ\text{C}$ . The broad peak at the extreme left shows the thermal conductivity cell response during pyrolysis; the subsequent sharp peak shows the response during flash release of the products from the pyrolyzer trap.

Figure 5. Schematic of the coupling of the gas chromatograph with the fast scan vapor-phase infrared spectrophotometer and the computer.

Figure 6. Mass chromatogram of the volatile pyrolysis products of low density polyethylene ( $\bar{M}_w = 140,000$ ,  $\bar{M}_n = 16,000$ ). Peak attenuations were X8, except for peaks 3 and 4 in the Freon-115 channel for which the attenuations were X64. Columns (SE-30) were program-heated at 5°C/min from 30-300°C.

Figure 7. Mass chromatogram of the volatile pyrolysis products of isotactic polypropylene (~100% isotactic,  $\bar{M}_w = 284,000$  and  $\bar{M}_w/\bar{M}_n = 2.98$ ). In the CO<sub>2</sub> channel peak attenuations were X32 for peak 9, and X8 for all others. In the Freon-115 channel attenuations were X128 (Peak 1), X64 (peaks 2 through 5), X32 (peaks 6 through 11), and X8 thereafter. Columns (SE-30) were program-heated at 7°C/min from 30° to 300°C.

Figure 8. Mass chromatogram of the volatile pyrolysis products of polyisobutylene ( $\bar{M}_v = 120,000$ ). Peak attenuations were X8 in both channels, except for peaks 1 and 2 for which the setting was X128. Columns (SE-30) were program-heated at 5°C/min from 30° to 300°C.

Figure 9. Comparison of the CO<sub>2</sub> channels of the mass chromatograms of volatile pyrolysis products of poly(butene-1-sulfone) (top), poly(pentene-1-sulfone) (middle), and poly(hexene-1-sulfone) (bottom).

Figure 10. Thermal conductivity cell responses during pyrolysis of atactic ( $\bar{M}_w = 105,000$ ,  $\bar{M}_w/\bar{M}_n = 2.15$ ), syndiotactic (70% syndio,  $\bar{M}_w = 83200$ ,  $\bar{M}_w/\bar{M}_n = 1.33$ ) and isotactic (91.5% iso,  $\bar{M}_w = 2,780,500$ ,  $\bar{M}_w/\bar{M}_n = 2.28$ ) polymethylmethacrylates.

Figure 11. Mass chromatogram of the volatile products of pyrolysis formed in the temperature range 200-500°C of atactic (A), syndiotactic (B) and isotactic (C) polymethylmethacrylates. Only the CO<sub>2</sub> channel outputs are shown for the syndiotactic and isotactic polymers. Columns (Dexsil-300) were program-heated at 5°C/min from 30 to 300°C. Peak attenuations were X64 for peak 1, and X8 for all others. The figure at the bottom (D) shows the CO<sub>2</sub> channel output for the products formed from the isotactic polymer in the temperature range 420 to 500°C.

Figure 12. Thermal conductivity cell responses during pyrolysis of atactic ( $\bar{M}_w = 359,950$ ,  $\bar{M}_w/\bar{M}_n = 2.03$ ), syndiotactic ( $\bar{M}_w = 20940$ ,  $\bar{M}_w/\bar{M}_n = 1.49$ ) and isotactic (82% iso,  $\bar{M}_w = 300,000$ ,  $\bar{M}_w/\bar{M}_n = 1.13$ ) poly(tertiarybutylmethacrylate)s.

Figure 13. Mass chromatogram of the volatile products of pyrolysis of isotactic poly(tertiarybutylmethacrylate). Columns (Dexsil-300) were program-heated at 5°C/min from 30° to 300°C. Peak attenuations were X8, except for peak 1 for which the settings were X128 in the Freon channel and X16 in the CO<sub>2</sub> channel.

Figure 14. Mass chromatogram of the volatile products of pyrolysis of anionic polystyrene ( $\bar{M}_n = 2100$ ,  $\bar{M}_w/\bar{M}_n < 1.1$ ) with n-butyl ends. Columns (Dexil-300) were program-heated at 5°C/min from 30° to 300°C.



Figure 15. Mass chromatogram of the volatile products of pyrolysis of thermal polystyrene ( $\bar{M}_n = 1880$ ,  $\bar{M}_w/\bar{M}_n = 1.01$ ). Columns (Dexsil 300) were program-heated at 5°C/min from 30° to 300°C.

Figure 16. Mass chromatogram of the volatile products of pyrolysis of anionic polystyrene ( $\bar{M}_n = 10500$ ,  $\bar{M}_w/\bar{M}_n = 1.15$ ) with  $\alpha$ -naphthyl ends. Columns (Dexsil-300) were program-heated at 5°C/min from 30° to 300°C.

Figure 17. Thermal history before and during pyrolysis of 1,4-polybutadienes. Solid lines for sample with 87% trans,  $\bar{M}_w = 67300$ ,  $\bar{M}_w/\bar{M}_n = 2.18$ ; and dashed lines for sample with 0% trans,  $\bar{M}_w = 633,500$ ,  $\bar{M}_w/\bar{M}_n = 3.62$ .

Figure 18. Mass chromatogram of the volatile pyrolysis products of high cis 1,4-polybutadiene (0% trans,  $\bar{M}_w = 633,500$ ,  $\bar{M}_w/\bar{M}_n = 3.62$ ) formed during 15 percent weight loss. Columns (Dexsil-300) were program-heated at 4°C/min from 30° to 300°C. Peak attenuations for the Freon-115 channel were X64 for peak 1, and X16 thereafter. That for the CO<sub>2</sub> channel was X8 for all peaks.

Figure 19. Mass chromatogram of the volatile pyrolysis products of high trans 1,4-polybutadiene (87% trans,  $\bar{M}_w = 67300$ ,  $\bar{M}_w/\bar{M}_n = 2.18$ ) formed during 15 percent weight loss. Columns (Dexsil-300) were program-heated at 4°C/min from 30° to 300°C. Peak attenuations for the Freon-115 channel were X64 for peak 1, and X16 thereafter. That for the CO<sub>2</sub> channel was X8 for all peaks.

Figure 20. "On-the-fly" infrared spectra of peaks 1,2,5 and 9 of the pyrogram shown in Figure 19.

Figure 21. Volatile degradation products from 1,4-polybutadienes: amounts of main products versus trans content (15 percent weight loss).

- 1,3-butadiene
- ① cyclopentene
- ① 1,3 cyclohexadiene
- 4-vinyl-1-cyclohexene
- △ cyclopentene + 4-vinyl-1-cyclohexene.

[Numbers identify samples (10)].

Figure 22. Mass chromatogram of the volatile pyrolysis products of high trans 1,4-polybutadiene (87% trans,  $\bar{M}_w = 67300$ ,  $\bar{M}_w/\bar{M}_n = 2.18$ ) formed during 15 percent to 100 percent weight loss. Columns (Dexsil-300) were program-heated at 4°C/min from 30° to 300°C. Peak attenuations for the Freon-115 channel were X64 for peak 1, and X16 thereafter. That for the CO<sub>2</sub> channel was X8 for all peaks.




Fig. 1








Figure 2

Figure 3





**Figure 4**






Fig. 5




Fig. 6




Fig. 7






Fig. 8




Fig. 9




Fig. 10






Fig. 11




Fig. 12




Fig. 13






Fig. 14




Fig. 15





Fig. 16








Figure 18




Figure 19






Figure 20




Figure 21




Figure 22



REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM								
1. REPORT NUMBER Technical Report No. 16	2. GOVT ACCESSION NO. (14) TR-26	3. ASSIGNMENT'S CATALOG NUMBER								
4. TITLE (and Subtitle) Pyrolysis-Molecular Weight Chromatography-Vapor Phase Infrared Spectrophotometry: An On-Line System for Analysis of Polymers	5. TYPE OF REPORT & PERIOD COVERED Technical Report, April 1978-April 1979	6. PERFORMING ORG. REPORT NUMBER								
7. AUTHOR(s) E. Kiran and J. K. Gillham	8. CONTRACT OR GRANT NUMBER(s) (15) N00014-76-C-0200									
9. PERFORMING ORGANIZATION NAME AND ADDRESS Polymer Materials Program Department of Chemical Engineering Princeton University, Princeton, NJ 08544	10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS Task No. NR 356-504									
11. CONTROLLING OFFICE NAME AND ADDRESS Office of Naval Research 800 North Quincy St. Arlington, VA 22217	12. REPORT DATE (11) May 1979	13. NUMBER OF PAGES 51								
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office) (12) 57 p.	15. SECURITY CLASS. (of this report)	15a. DECLASSIFICATION/DOWNGRADING SCHEDULE								
16. DISTRIBUTION STATEMENT (of this Report) Approved for Public Release; Distribution Unlimited										
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report) (10) Erdogan/Kiran John K./Gillham										
18. SUPPLEMENTARY NOTES										
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) <table border="0"> <tr> <td>Pyrolysis</td> <td>Poly olefin sulfones</td> </tr> <tr> <td>Mass Chromatography</td> <td>Poly methacrylates</td> </tr> <tr> <td>Vapor Phase Infrared Spectrophotometry</td> <td>Polystyrenes</td> </tr> <tr> <td>Polyolefins</td> <td>Polybutadienes</td> </tr> </table>			Pyrolysis	Poly olefin sulfones	Mass Chromatography	Poly methacrylates	Vapor Phase Infrared Spectrophotometry	Polystyrenes	Polyolefins	Polybutadienes
Pyrolysis	Poly olefin sulfones									
Mass Chromatography	Poly methacrylates									
Vapor Phase Infrared Spectrophotometry	Polystyrenes									
Polyolefins	Polybutadienes									
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) <p>An on-line pyrolysis system has been developed which comprises a programmable pyrolyzer, a thermal conductivity cell, a trap, a mass chromatograph, a conventional gas chromatograph, a vapor-phase infrared spectrophotometer, and a digital computer. The system has proved powerful for analysis of polymers as has been shown by results obtained with polyolefins, polyolefin-sulfones, polymethacrylates, polystyrenes and polybutadienes. This report provides a brief description of the system and reviews results obtained with these polymers.</p>										

DD FORM 1 JAN 73 1473

EDITION OF 1 NOV 65 IS OBSOLETE  
S/N 0102-LF-014-6601

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

404 377

slv

TECHNICAL REPORT DISTRIBUTION LIST, GEN

	<u>No.</u> <u>Copies</u>		<u>No.</u> <u>Copies</u>
Office of Naval Research 800 North Quincy Street Arlington, Virginia 22217 Attn: Code 472	2	Defense Documentation Center Building 5, Cameron Station Alexandria, Virginia 22314	12
ONR Branch Office 536 S. Clark Street Chicago, Illinois 60605 Attn: Dr. George Sandoz	1	U.S. Army Research Office P.O. Box 1211 Research Triangle Park, N.C. 27709 Attn: CRD-AA-IP	1
ONR Branch Office 715 Broadway New York, New York 10003 Attn: Scientific Dept.	1	Naval Ocean Systems Center San Diego, California 92152 Attn: Mr. Joe McCartney	1
ONR Branch Office 1030 East Green Street Pasadena, California 91106 Attn: Dr. R. J. Marcus	1	Naval Weapons Center China Lake, California 93555 Attn: Dr. A. B. Amster Chemistry Division	1
ONR Area Office One Hallidie Plaza, Suite 601 San Francisco, California 94102 Attn: Dr. P. A. Miller	1	Naval Civil Engineering Laboratory Port Hueneme, California 93401 Attn: Dr. R. W. Drisko	1
ONR Branch Office Building 114, Section D 666 Summer Street Boston, Massachusetts 02210 Attn: Dr. L. H. Peebles	1	Professor K. E. Woehler Department of Physics & Chemistry Naval Postgraduate School Monterey, California 93940	1
Director, Naval Research Laboratory Washington, D.C. 20390 Attn: Code 6100	1	Dr. A. L. Slafkosky Scientific Advisor Commandant of the Marine Corps (Code RD-1) Washington, D.C. 20380	1
The Assistant Secretary of the Navy (R,E&S) Department of the Navy Room 4E736, Pentagon Washington, D.C. 20350	1	Office of Naval Research 800 N. Quincy Street Arlington, Virginia 22217 Attn: Dr. Richard S. Miller	1
Commander, Naval Air Systems Command Department of the Navy Washington, D.C. 20360 Attn: Code 310C (H. Rosenwasser)	1	Naval Ship Research and Development Center Annapolis, Maryland 21401 Attn: Dr. G. Bosmajian Applied Chemistry Division	1
		Naval Ocean Systems Center San Diego, California 91232 Attn: Dr. S. Yamamoto, Marine Sciences Division	1

Encl 1

TECHNICAL REPORT DISTRIBUTION LIST, 356A

	<u>No.</u> <u>Copies</u>		<u>No.</u> <u>Copies</u>
Dr. Stephen H. Carr Department of Materials Science Northwestern University Evanston, Illinois 60201	1	Picatinny Arsenal SMUPA-FR-M-D Dover, New Jersey 07801 Attn: A. M. Anzalone Building 3401	1
Dr. M. Broadhurst Bulk Properties Section National Bureau of Standards U.S. Department of Commerce Washington, D.C. 20234	2	Dr. J. K. Gillham Princeton University Department of Chemistry Princeton, New Jersey 08540	1
Dr. T. A. Litovitz Department of Physics Catholic University of America Washington, D.C. 20017	1	Douglas Aircraft Co. 3855 Lakewood Boulevard Long Beach, California 90846 Attn: Technical Library CI 290/36-84 AUTO-Sutton	1
Dr. R. V. Subramanian Washington State University Department of Materials Science Pullman, Washington 99163	1	Dr. E. Baer Department of Macromolecular Science Case Western Reserve University Cleveland, Ohio 44106	1
Dr. M. Shen Department of Chemical Engineering University of California Berkeley, California 94720	1	Dr. K. D. Pae Department of Mechanics and Materials Science Rutgers University New Brunswick, New Jersey 08903	1
Dr. V. Stannett Department of Chemical Engineering North Carolina State University Raleigh, North Carolina 27607	1	NASA-Lewis Research Center 21000 Brookpark Road Cleveland, Ohio 44135 Attn: Dr. T. T. Serofini, MS-49-1	1
Dr. D. R. Uhlmann Department of Metallurgy and Material Science Center for Materials Science and Engineering Massachusetts Institute of Technology Cambridge, Massachusetts 02139	1	Dr. Charles H. Sherman, Code TD 121 Naval Underwater Systems Center New London, Connecticut	1
Naval Surface Weapons Center White Oak Silver Spring, Maryland 20910 Attn: Dr. J. M. Augl Dr. B. Hartman	1	Dr. William Risen Department of Chemistry Brown University Providence, Rhode Island 02192	1
Dr. G. Goodman Globe Union Incorporated 5757 North Green Bay Avenue Milwaukee, Wisconsin 53201	1	Dr. Alan Gent Department of Physics University of Akron Akron, Ohio 44304	1



TECHNICAL REPORT DISTRIBUTION LIST, 356A

	<u>No.</u> <u>Copies</u>		<u>No.</u> <u>Copies</u>
Mr. Robert W. Jones Advanced Projects Manager Hughes Aircraft Company Mail Station D 132 Culver City, California 90230	1	Dr. T. J. Reinhart, Jr., Chief Composite and Fibrous Materials Branch Nonmetallic Materials Division Department of the Air Force Air Force Materials Laboratory (AFSC) Wright-Patterson Air Force Base, Ohio 45431	1
Dr. C. Giori IIT Research Institute 10 West 35 Street Chicago, Illinois 60616	1	Dr. J. Lando Department of Macromolecular Science Case Western Reserve University Cleveland, Ohio 44106	1
Dr. M. Litt Department of Macromolecular Science Case Western Reserve University Cleveland, Ohio 44106	1	Dr. J. White Chemical and Metallurgical Engineering University of Tennessee Knoxville, Tennessee 37916	1
Dr. R. S. Roe Department of of Materials Science and Metallurgical Engineering University of Cincinnati Cincinnati, Ohio 45221	1	Dr. J. A. Manson Materials Research Center Lehigh University Bethlehem, Pennsylvania 18015	1
Dr. L. E. Smith U.S. Department of Commerce National Bureau of Standards Stability and Standards Washington, D.C. 20234	1	Dr. R. F. Helmreich Contract RD&E Dow Chemical Co. Midland, Michigan 48640	1
Dr. Robert E. Cohen Chemical Engineering Department Massachusetts Institute of Technology Cambridge, Massachusetts 02139	1	Dr. R. S. Porter University of Massachusetts Department of Polymer Science and Engineering Amherst, Massachusetts 01002	1
Dr. David Roylance Department of Materials Science and Engineering Massachusetts Institute of Technology Cambridge, Massachusetts 02039	1	Professor Garth Wilkes Department of Chemical Engineering Virginia Polytechnic Institute and State University Blacksburg, Virginia 24061	1
Dr. T. P. Conlon, Jr., Code 3622 Sandia Laboratories Sandia Corporation Albuquerque, New Mexico	1	Dr. Kurt Baum Fluorochem Inc. 6233 North Irwindale Avenue Azusa, California 91702	1
Dr. Martin Kaufmann, Head Materials Research Branch, Code 4542 Naval Weapons Center China Lake, California 93555	1	Professor C. S. Paik Sung Department of Materials Sciences and Engineering Massachusetts Institute of Technology Cambridge, Massachusetts 02139	1